



Ready-to-Use

General Virucide, Bactericide, Tuberculocide, Fungicide, Sanitizer

One-step cleaner disinfectant formulated with patented hydrogen peroxide technology to deliver faster, more effective cleaning performance.

Meets OSHA's bloodborne pathogen standards for HBV and HIV.

Broad Spectrum Germicidal Action

BACTERICIDAL and *VIRUCIDAL in one minute. Kills HBV, HCV, HIV-1 (AIDS virus) in one minute. TUBERCULOCIDAL in five minutes. Kills fungi in ten minutes.

Improved Cleaning

The patented technology** contains a high performance surfactant system that delivers faster and better cleaning results. No rinsing is required in most applications for improved productivity.

Environmentally Responsible

No VOC's (Volatile Organic Compounds), no added fragrance, no NPE's (nonylphenol ethoxylates). The active ingredient breaks down to water and oxygen.

Non-Corrosive Formula

Can be used to clean most hard, non-porous environmental surfaces with confidence. User Safety—non-irritating to skin and eyes and contains no substances which at the given concentration, are considered to be hazardous to health per OSHA standards.



Johnson wax
PROFESSIONAL



OXIVIR™/MC Tb

General Virucide, Bactericide, Tuberculocide, Fungicide, Sanitizer

Profile

WHAT IS IT?	A revolutionary, patented one-step cleaner disinfectant for spray and wipe cleaning applications.
WHAT DO I USE IT ON?	Use daily on the following hard, non-porous environmental surfaces: vinyl, plastic surfaces, glazed ceramic, glazed porcelain, chrome, laminated surfaces (and baked enamel surfaces associated with floors), walls, ceilings, tables, chairs, countertops, telephones, fixtures, glazed tile, toilet exteriors, urinal exteriors, sinks found in health care facilities, hospitals, patient rooms, operating rooms, shower rooms, and locker rooms. It can be used to pre-clean and disinfect hospital items: wheelchairs, hospital patient bedrails and linings, patient bed linings, vinyl mattress covers, wash basins, bed pans, non-critical medical equipment surfaces – almost any washable food and non-food contact surface where disinfection is required. A potable water rinse is required when disinfecting food contact surfaces. Do not use on glasses, dishes and utensils.
HOW DO I USE IT?	Pre-clean heavily soiled areas. Apply solution by spray to hard, non-porous environmental surfaces. All surfaces must remain wet for 1 minute to kill bacteria, HIV-1, HBV and HCV. Use a 5 minute contact time for Tb and a 10 minute contact time for fungi. Wipe surfaces dry, rinse, or allow to air dry.
NOTE:	Rinsing is not necessary unless floors are to be coated with finish or restorer. All food contact surfaces must be rinsed with potable water. Do not use on glassware, utensils, or dishes.

Specifications

Active Ingredients: Hydrogen Peroxide 0.5%; **EPA Reg. No.:** 70627-56;
Color / Form: Clear colorless liquid; **pH:** 2.5 – 3.5; **Fragrance:** No fragrance
Cleaning Performance: Excellent; **Phosphates:** None; **Specific Gravity:** (water = 1) 1.01

Efficacy Summary

This product effectively kills the following on hard, non-porous environmental surfaces in the presence of 5% serum load and 1 minute contact time (unless otherwise noted)

Bacteria:

Pseudomonas aeruginosa, Escherichia coli 0157:H7, Acinetobacter baumannii, Staphylococcus aureus, Staphylococcus aureus MRSA, Salmonella choleraesuis, Enterococcus faecalis VRE

*Viruses:

Poliovirus Type 1, Strain Brunhilde
Human Coronavirus
Herpes Simplex Virus, Type 2 (HSV-2)
Hepatitis C Virus (HCV)
HIV-1 (AIDS Virus), Strain HTLV-IIIIB (HIV-1)
Influenza A/Hong Kong, Rhinovirus Type 37, Strain 151-1, Hepatitis B Virus (HBV)
Feline Calicivirus, Strain F9 Herpes Simplex Virus, Type 1 (HSV-1) Rotavirus WA

TUBERCULOCIDAL Activity: (5 minute contact time)

Mycobacterium bovis (Tb) (OT 451C150)

FUNGI – (10 minute contact time)

Tricophyton mentagrophytes

Non-Food Contact Surface Sanitizer – (30 second contact time)

Klebsiella pneumoniae, Staphylococcus aureus, Escherichia coli 0157:H7, Salmonella choleraesuis, Staphylococcus aureus MRSA, Pseudomonas aeruginosa, Enterococcus faecalis VRE

Safety Reminder

Please make sure your employees read and understand the product label and Material Safety Data Sheet before using this product. The label contains directions for use, and both the label and MSDS contain hazard warnings, precautionary statements and first aid procedures. MSDS's are available online at www.johnsondiversey.com or by calling 1-888-352-2249.

HMIS/NFPA HAZARD RATING

Health 0
Flammability 0
Reactivity 0



Packaging

PACK / SIZE	PRODUCT CODE
12 / 32 fl. oz. containers	4277285

Johnson wax
PROFESSIONAL

FOR COMMERCIAL & INDUSTRIAL USE ONLY

For further information on special applications or other Johnson Wax Professional products, contact your JohnsonDiversey representative, call us at 1-800-558-2332 or visit our web-site at www.johnsondiversey.com.